


ACURIS™ Life



Open your ears to
a Life of better hearing.

Siemens Hearing Instruments
Better hearing is better living.

SIEMENS



Introducing

ACURIS Life

Open ear hearing instruments have enabled many people with mild to moderate hearing loss to enjoy fuller lives. They're easy to wear. Easy to use. And barely noticeable. The only thing they lack is knockout performance. That is, until now. Introducing ACURIS™ Life. It's a quantum leap forward in open ear hearing instrument technology. With a level of performance no other open solution can match.



ACURIS Life

Nothing else comes close.



The tiny ACURIS Life instrument is only 1 7/16 inches long (30 mm) and weighs less than 1/10 of an ounce (2.7 grams). When placed behind the ear it's virtually invisible.

(Shown actual size)

A subtle look.

ACURIS Life's open ear solution utilizes the comfortably soft, flexible, and skin-tunesure™ silicone dome to help reduce the bothersome occlusion effect — the echo you may feel when wearing a hearing aid. With ACURIS Life, no one will notice the hearing aid is there.

Easy to use.

ACURIS Life's instrument couldn't be simpler. No buttons, no dials, no instrument to fiddle with. Instead, you can use the optional Pocket remote. This small, sleek remote lets you adjust the volume and power of your ACURIS Life instrument.



ACURIS Life is the only open ear instrument with the optional Pocket remote, which lets you have complete control conveniently at your fingertips.

Unmatched performance.

Available ACURIS Life instruments come with features you won't find in any other open system, including Advanced High Speed Feedback Cancellation and Advanced Digital Signal Processing with sophisticated Digital Noise Suppression and Digital Speech Enhancement.

Should you desire to wear two instruments, only ACURIS Life offers e2e wireless, the engineering breakthrough that enables both instruments to communicate with each other and function as one, holistic hearing system that adjusts automatically for optimal hearing in any environment. Whether you wear ACURIS Life in one or both ears, you'll benefit from hearing performance unmatched by any other open instrument.

To learn more about ACURIS Life, speak with your Hearing Care Professional.

ACURIS Life is available at:

Siemens Hearing Instruments, Inc.
P.O. Box 1397
Piscataway, NJ 08855-1397
www.usa.siemens.com/hearing

Siemens Hearing Instruments
A Division of Siemens Canada Limited
320 Pinebush Road, Cambridge, Ontario N3C 2V3

Information in this brochure is subject to change without notice. Copyright © 2005 Siemens Hearing Instruments, Inc. All rights reserved.

Hearing instruments may help people hear better, but no hearing instrument can solve every problem or restore hearing back to normal.